

SPOT PathTracker™

Bulk Barcode Reader

PathTracker is capable of reading 150+ linear/2D barcodes in under 25 seconds, saving your lab time and money. It provides advanced sample tracking and loss prevention. PathTracker LIS integration transfers scanned results directly to your LIS.



Key Values:

- Scan All Transfers of:
 - Cassettes, Slides and Vials
- Use Existing Standard Baskets
- Easy to Use Graphic Interface
- Stand-Alone or LIS Linked
- SQL Data Base Archives of Scan Data

PathTracker Features:

- Scan 150 Cassettes in <25s
- Scan 20 slides in <15s
- Automatic Manifest Printout
- Receiving Reconciliation
- Late Sample Auto-Notification
- Duplicate Sample Flag
- Stat Sample Flagging
- Correction Tools for Damaged Barcodes
- Sample Tracking Report




PathSuite™

Integrated Imaging Environment

The Vision of the Future: Connected Solutions

Building a vision of the future is essential to ensuring that each step supports the final integrated solution. PathSuite™ provides an over arching vision that supports streamlined workflows for the task at hand built into a connected environment that supports collaboration.



SPOT PathSuite™ Portfolio

See all of our Pathology Solutions at:
spotimaging.com