



PathSuite™ Software


Imaging Software Optimized for Pathology




PathSuite™ | Easy to Learn - Easy to Use

PathSuite provides the user with a simple icon driven user interface with one click access to essential controls. This makes PathSuite intuitive to learn and a pleasure to use. PathSuite supports both gross and microscopy imaging with a common user interface.

Gross Imaging



Microscopy Imaging



It's easy!

PathSuite™ | Streamlined Workflow

PathSuite's streamlined workflow saves time on every image through the tight integration of instrument hardware coupled with the automation of Login, Preview and Archiving functions. The results add up to hours of savings per day.

Log In

PathSuite provides quick, accurate validation of your image with case Log in.



File Server



LIS



Cloud

Archive

A single click captures and saves images to multiple archive locations, adding annotations as needed.




Preview

PathSuite makes image composition easy with live preview of your camera adjustments.

Get your work done!


PathSuite™ | Measure and Annotate

PathSuite provides a full suite of annotation and measurement tools for detailing your images without transferring images between programs or computers.




Calibration Mark

Calibration marks add a dimensional reference to your image.




Measurement

PathSuite's pre-calibrated images allow for easy measurement of length, area and angle.




Standard Annotations

A full set of text and drawing tools ensure you have the arrow text or oval to quickly and fully detail your images for consultation, publication or standard reporting.




Extended Annotations

Speed repetitive pathology annotations with custom text imprints, as well as auto-numbered labels, sections and blocks. Additionally, auto-imprint standard annotations with each image archived.



Style

Annotation style preferences are easily set or changed with control over numbering protocol, font, size, color, fill, etc...



Case Management Tools

Easily select, duplicate and move images to support your "after hours" and case management work.


Express Yourself

PathSuite™ | Telepathology

PathSuite provides simple, secure connections and fully interactive live consults in a practical and affordable format. Now every workstation and computer can be virtually connected making telepathology a routine part of your teams work practices. Amazingly, this increase in effectiveness comes with an increased efficiency from the elimination of unnecessary travel.

Save Time


PathSuite telepathology solutions save you time on each trip you make whether it be down the hall, down the street or across the world. It eliminates the time and frustration associated with the trip and the unaccounted for impromptu interactions that occur with each trip. Telepathology lets you connect, consult and return to work without leaving your desk, it's that simple.



(*Includes: Scheduling time, hospital walk time, impromptu discussion time, parking time, congested traffic time, gas stop time, early arrival buffer time, check in security....]

Connect

Colleague meetings are initiated by selecting the team member from the connection menu. The colleague receives an invitation prompt on their screen and joins the meeting with a simple click. Non-colleague attendees are easily added to meetings via e-mail invites. Once accepted, meetings are secured using guest passwords and data encryption, supporting HIPPA requirements.




Consult

PathSuite ensures a rich collaborative meeting allowing colleagues to easily view images, discuss interpretations and add annotations. The shared cursor and annotation tools allow each attendee to direct the discussion as needed.



Document

Discussion images are easily saved to the case with a click of a button preserving results for later reference and follow up.



Be Effective, Be Efficient

Return on Your Investment

Our install base has garnered dramatic time savings on image acquisition across multiple pathology work tasks resulting in impressively short ROI periods, followed by ongoing savings year after year.

SAVE TIME!

Routine Pathology Slide Examination

Time Saved per Slide

Slides per Day

Time Saved per Day

Time Saved per Year

2 Min

X

60

2 Hrs

or

520 Hrs

Routine Gross Dissection

Time Saved per Specimen

X

Specimens per Day

Time Saved per Day

Time Saved per Year

10 Min

12

$$=$$

2 Hrs

or

520 Hrs

Increase Consultations



Increasing intraoperative and needle biopsy adequacy consultations provides savings for the organization and the physician by eliminating unreimbursed secondary procedures and increasing patient satisfaction. Pathologists and Pathologist Assistants benefit from the increased efficiency of providing their consults virtually, avoiding the unproductive trip otherwise required.

[illegible]

PathSuite™ | The Smart Choice

Professional Documentation

Consistent documentation of your work provides the foundation for quality assurance and process improvement. It is also the basis of how you are perceived by your patients, referring physicians and by the legal establishment.



Easy to Implement

SPOT Imaging Solutions has developed turnkey solutions with configuration options to support your facility's unique hardware and software infrastructure. The result is an optimized work environment that is easily implemented into your existing infrastructure.

"The installation was quick and we were imaging the next day!"

*Satisfied Pathology Assistant,
Canada*

Easy to Support

You are never alone with PathSuite. InSpect™ technical service. Our 47 year tradition of standing behind our equipment ensures your uptime performance with on-line, on-phone or on-site assistance.

"We had a different system that was sitting in a box because we couldn't use it, PathSuite is so easy to use... Thanks SPOT Imaging!"

Medical Center, MI

Pathology Solutions to Support Your Team

Microscope Imaging

- Routine Anatomical Pathology
- Frozen Section Consults
- Needle Biopsy Consults
- Tumor Board
- Teaching and Training

Macro Imaging

- Routine Gross Dissection
- Frozen Section Consults
- Teaching and Training

PathSuite ...The Standard Imaging Interface for the Pathology Laboratory

PathSuite™ | Industry Leading Performance

Administrators...

"...the savings we achieve with the PathSuite software makes it a no brainer investment justification. We feel that we are now a state of the art pathology operation..."

University, NC

User's...

"So intuitive and easy to use! ... now that we have SPOT, we get our work done quickly and without frustration! It works great! "

University, GA

IT/Technical Support...

"...SPOT stands behind their products with on line support, repair services and reliable products ... They've been there for us! "

Medical Center, PA

Some of our Customer Locations...

University of Miami	Carillion Canada	Tricare Reference Laboratories
West Virginia University Hospital	Central Hospitals of Montreal	Temple University Hospital
University of Georgia	Cedars-Sinai Medical Center	Duke University
Laval University-Eastern Quebec Telepathology	Maine Medical Center	Kaiser Hospital
Broward Health	The University of Florida	Nycom Inc.
California Pacific Medical Center	Reading Hospital	Virginia Hospital Center
Orlando Regional Health Systems	Defense Finance & Acct. Service/DFAS Rome	University of Pittsburgh Medical Center
Texas Children's Hospital	Children's Healthcare of Atlanta	University of Pennsylvania Health System
Engineering Systems Inc.	Presbyterian Healthcare Services of NM	Alberta Provincial Healthcare
Michigan State University	Appalachian Reg Healthcare System	
Children's Hospital LA		

Call SPOT Imaging Solutions for experienced advice on how to dramatically improve your pathology lab with SPOT imaging and telepathology



6540 Burroughs Avenue • Sterling Heights • Michigan 48314-2133 USA
Toll Free USA: **1-866-604-SPOT** • Phone: **+1 586-731-6000** • Fax: **+1 586-731-6469**
Web: **www.spotimaging.com**