

SPOTBroadcast: Real Time Telepathology Broadcasting

SPOTBroadcast captures video in high definition and provides cost effective, one-way streaming of images over the internet without a computer.

SPOTBroadcast enables gross specimens, autopsy, gram stains, and fine needle aspirations to be safely and securely broadcasted over the internet in real time to remote locations anywhere. Exceptional diagnostic quality HD images of gross specimens and glass slides can be shared conveniently for review and diagnosis.

SPOTBroadcast captures audio and video sources and streams the input so that it can be shared remotely at full HD resolution. SPOTBroadcast can be used anywhere - it's portable with no software required.



SPOTBroadcast™




Share Images With Ease



Barcode



Footswitch



Speakerphone

Features	Benefits
High Resolution Camera	Diagnostic quality images with accurate color and fast frame rates are ideal for documentation, pathology reports, and live consultations.
Speakerphone	Dictate description and diagnoses
Hardware Broadcaster	Streams the HD video and audio to remote colleagues.
Barcode Scanner	Scanning in case information saves time and reduces errors.
Footswitch Control	Control the system without looking up from your microscope.
Portable	SPOTBroadcast is easy to move where diagnoses are needed.
Flexible Configuration	SPOTBroadcast is configurable to fit your specific needs and budget



SPOT Imaging Solutions™

A Division of Diagnostic Instruments, Inc.