



MagMinder™* | Practical

*Patent Pending


MagMinder™ : track your magnification settings




Key Features

- Fits Existing Microscopes
- Automatic Reading of Magnification
- Automatic Setting of Image Calibration
- Automatic Recording of Magnification to Image Memo
- Automatic Imprint of Calibration Mark


Fits Nearly Any Microscope




Magnification Recording



Auto-Calibrated Measurements



Auto-Imprint of Calibration Mark



MagMinder™* | The Smart Choice

*Patent Pending


MagMinder automatically tracks, records and imprints your images with the current magnification, calibration setting and reference mark. These automatic functions can save up to 2 minutes per slide over manual setting, recording and imprinting. This can have a significant impact in production environments.



Time saved per slide		Slides per day		Time saved per day		Time saved per year
2 Min.	X	60	=	2 Hrs	or	520 Hrs

The Time is Yours ...You Decide How to Use It!

Recording Assurance




MagMinder tracks your magnification settings so you can concentrate on your research. When needed, your image acquisition data is just a click away in the SPOT Image Information Dialog: date, exposure settings and now your magnification.

Upgrade and Save



Upgrading your microscope with SPOT automation accessories saves thousands of dollars over the purchase of a replacement microscope. SPOT accessories modular designs ensure quick, easy service and maintenance.


Easy to Install



MagMinder installs in minutes on nearly any microscope and runs conveniently under SPOT Advanced Software with a click of the acquire button. Installation can be performed by a SPOT Imaging Specialist or by the customer.

Go...

Easy to Learn



Once installed, operating MagMinder is as simple as clicking the acquisition button. The MagMinder learns your magnifications in its first session via prompts and then runs seamlessly from then on. The simple settings dialog makes calibration imprints and magnification resetting a snap. MagMinder makes it simple.



6540 Burroughs Avenue
Sterling Heights • MI • 48314-2133
phone: 586.731.6000 • fax: 586.731.6469
website: www.spotimaging.com

SPOT MagMinder

Call SPOT Imaging Solutions for expert advice on how you can transform your pathology environment with digital imaging and telepathology