

Need to show your students detailed dissections without enough room or time?

Make It
EASY
to See!

PathMobile™ easily positions over your dissection, holds the camera for hours and projects wirelessly to your Lab's monitors!

PathMobile™ Makes Teaching Anatomy Easy...

- Safe and Easy Positioning of Camera
- 20 Mp Image Capture and Auto Archive
- Wireless 720p HD Live Video Projection to Room Monitor
- 720p HD Audio-Video Recording
- Internet Broadcasting of Your Lecture
- Full Body to 2.5 cm FOV
- Easy to use Annotation Tools
- Add Archived Reference Images to Supplement Your Lecture
- 6 Hour Cordless Battery Operation


Now All your students can see your dissections in the lab, across campus and beyond...



PathMobile™


"Roll Around" Imaging Workstation

In the Lab




Project Video Wirelessly to Room Monitors ensuring every student has a clear view.

Remote Locations




Use **PathCast** Internet Broadcasting to support locations outside of the lab.

Record



Video record presentations and lectures for tutorials and reference of key dissections.


Complete Presentation



Easily transition between live video and captured images. Use digital zoom and annotations to make your point.




Max FOV



Capture the Details!

Max Magnification



Contact Us for a **Virtual** demo...



PathMobile's 6-axis positioning and 34 x zoom range ensure you and your students see the field of view required. From 205 cm whole body views to close up 2 cm magnified views.

6540 Burroughs Ave. • Sterling Heights, MI 48314
 Email: Info@spotimaging.com • Toll Free: +1 866-604-SPOT
 Direct: +1 586-731-6000
 PN: University Anatomy Flyer