



# PATHARCHIV™

Powered by Lab Improvements

PATH ARCHIV

## The Next Step in File Room Management

### Automated digital block and slide storage

- Eliminates sample losses from typos and misfiling
- Frees-up staff time required to manage file room
- Tracks samples in archives and during loan out
- 2.3 times the slide density of competing systems
- Uses standard-size storage format containers



**PathArchiv-B**  
Block Filing



**PathArchiv-S**  
Slide Filing

**\$140K**  
Per Year Savings!  
287,000 Blocks  
575,000 Slides



**PathArchiv** automates the recording, filing and loans of blocks and slides, dramatically reducing lost samples as well as the staff required to manage the file room.





# PathArchiv Specifications

	PathArchiv-S	PathArchiv-B
Computer OS	Win 10 Pro	Win 10 Pro
USB 3.0 Ports	2	2
LAN Port	1	1
Keyboard and Mouse	Standard Wireless	Standard Wireless
Handheld Barcode Reader	Optional	Optional
<b>Display</b>		
Touch Screen	Full color LCD	Full color LCD
<b>Power / Environment</b>		
Current	1.0 amp	1.0 amp
Power Supply	120-240 VAC	120-240 VAC
Frequency	50/60 Hz	50/60 Hz
Operating Temperature	41° - 86° F / +5 - +30°C	41° - 86° F / +5 - +30°C
Certifications	CE, FCC Class A, RoHS CSA/UL Electrical Certification	CE, FCC Class A, RoHS CSA/UL Electrical Certification
<b>Instrument Size</b>		
Height	17" H (431.8 mm)	17.5" H (444.5 mm)
Width	27" W (685.8 mm)	10.75" W (273.05 mm)
Depth	19" D (482.6 mm)	22.5" D (571.5 mm)
<b>Weight</b>		
Instrument	52 lb / 23.6 kg	34 lb / 15.4 kg
<b>Sample Size Capacity</b>		
Standard Slide Size	1.007" x 3.007" (25.6 mm x 76.4 mm)	
Standard Block Cassette Size		1.122" x 1.614" (28.5 mm x 41 mm)
<b>Sample Processing and Storage</b>		
Samples per Instrument	200	50
Samples per Magazine	200	50
Magazines per Storage Box	8	8
Samples per Storage Box	1,600	400
Storage Box Size	12" W (304.8 mm) 12" D (304.8 mm) 4" H (101.6 mm)	12" W (304.8 mm) 12" D (304.8 mm) 4" H (101.6 mm)
Magazine Size	2-5/8" W (66.675 mm) 5-5/8" D (142.875 mm) 3-1/8" H (79.375 mm)	2-5/8" W (66.675 mm) 11-1/8" D (282.575 mm) 1-1/4" H (31.75 mm)
<b>Sample Processing Rate (Samples per Hour)</b>		
	400 Slides/hr	1500 Blocks/hr


## PathArchiv-S Slide Filing




## PathArchiv-B Block Filing




### Slide storage magazine



### Block storage magazine



### Outer storage box



SPOT IMAGING™  
Take a closer look

586.731.6000 • toll free (US): 1-866-604-SPOT  
www.spotimaging.com • info@spotimaging.com